1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrume	nt serial N	0:	68-012337		Accuracy che	eck date/time:		08:21
			00 0 .200				00/00/2010	
			on check:		Interferen	t detected:		
	Second ca	llibrati	on check:		Interferen	t detected:		
Concenti	ration (%W	//V)	Solution lot#	ŧ	Soln to Simulator			
Number of Operator		ested	since last periodic de	termination of a	accuracy:			
				Diagnostics	check			
Prom F Ram F	Pass Pass Pass Jum 15351		Serial No. Pass Time clock Pass calibration Pass Download date/time:	Temp Pas Motor Pas Range Pas 06/05/2015	SS	Digital voltme	ter Pass	
	Reco	ord im	ported (Y/N) N	Sourc	e Diagnostic			
Instrume	nt serial N	0:	68-012337		Accuracy che	eck date/time:	06/05/2015	08:22
	First ca	libratio	on check: Not perform	ied	Interferen	t detected:	No	
	Second ca	llibrati	on check:Not perform	ied	Interferen	t detected:	No	
Concenti	ration (%W	//V)	Solution lot#	t:	Soln to Simulator			
Number of Operator		ested	since last periodic de	termination of a	accuracy: 030			
				Diagnostics	check			
Eeprom Prom Ram			Serial No. Time clock calibration	Temp Motor Range	Ir	nternal standa Prin Digital voltme	ter	
Check su	ım 50878 Reco		Download date/time: ported (Y/N) N	06/05/2015	15:40 e Calibration	_	d ID 8174	

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: 68-012337 Accuracy check date/time: 06/05/2015 08:28

First calibration check: .116 Interferent detected: No Second calibration check: .118 Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C: 33.99

Simulator Ser No. DR3707

Number of people tested since last periodic determination of accuracy: 030

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 40331 Download date/time: 06/05/2015 15:40 Record ID 8181

Record imported (Y/N) N Source Calibration

Instrument serial No: 68-012337 Accuracy check date/time: 06/09/2015 02:27

First calibration check: Interferent detected: Second calibration check: Interferent detected:

Concentration (%W/V) Solution lot#: Soln temp C:

Simulator Ser No.

2

Invalid sample

Number of people tested since last periodic determination of accuracy:

Operator

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 46858 Download date/time: 06/15/2015 09:51 Record ID 81G7

Record imported (Y/N) N Source Exception

This is a true, complete, and accurate report of the periodic determination of accuracy and diagnostic data from the specified Intoxilyzer 5000EN instrument obtained from the CMI COBRA database.

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No	o: 68-012337	A	Accuracy check date/time: 06	6/09/2015 03:06		
	llibration check: alibration check:		Interferent detected: Interferent detected:			
Concentration (%W	//V) Solution	lot#:	Soln temp C: Simulator Ser No.			
Invalid sample						
Number of people to	ested since last periodic	determination of ac	curacy:			
Operator						
		Diagnostics c	neck			
Eeprom	Serial No.	Temp	Internal standards	5		
Prom	Real Time clock	Motor	Printer			
Ram	Auto calibration	Range	Digital voltmeter			
Check sum 11989	Download date/ti		9:51 Record II	O 81H8		
Reco	ord imported (Y/N) N	Source	Exception			
Instrument serial No	o: 68-012337	A	Accuracy check date/time: 06	6/10/2015 03:14		
First ca	llibration check:		Interferent detected:			
Second ca	alibration check:		Interferent detected:			
Concentration (%W	//V) Solution	lot#:	Soln temp C: Simulator Ser No.			
Invalid sample						
Number of people to	ested since last periodic	determination of ac	curacy:			
Operator						
Diagnostics check						
Eeprom	Serial No.	Temp	Internal standards	5		
Prom	Real Time clock	Motor	Printer			
Ram	Auto calibration	Range	Digital voltmeter	r		
Check sum 44576	Download date/ti	mo: 06/15/2015 00	9:51 Record II	N 81I3		
	ord imported (Y/N) N		Exception	3 0 0		

This is a true, complete, and accurate report of the periodic determination of accuracy and diagnostic data from the specified Intoxilyzer 5000EN instrument obtained from the CMI COBRA database.

3

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial N	o: 68-0	012337		Accuracy	check date/time: 06/11/2015	11:33
First calibration check: Second calibration check:				Interferent detected: Interferent detected:		
Concentration (%W	//V)	Solution lot#:			n temp C: or Ser No.	
Invalid sample						
Number of people to	ested since	last periodic dete	rmination of	accuracy:		
Operator						
			Diagnostics	check		
Eeprom	Ser	ial No.	Temp		Internal standards	
Prom	Real Time		Motor		Printer	
Ram	Auto calib		Range		Digital voltmeter	
Check sum 13220	_	nload date/time: 0		09:51	Record ID 81J0	
Reco	ord importe	d (Y/N) N	Sour	ce Exception)	
Instrument serial N	o: 68-0	012337		Accuracy	check date/time: 06/13/2015	03:14
First ca	llibration ch	eck:		Interfer	ent detected:	
Second ca	alibration ch	ieck:		Interfer	ent detected:	
Concentration (%W	//V)	Solution lot#:			n temp C: or Ser No.	
Invalid sample						
Number of people to	ested since	last periodic dete	rmination of	accuracy:		
Operator						
Operator						
Operator		Г	Diagnostics	s check		
	Ser	Eial No.	Diagnostics Temp	s check	Internal standards	
Eeprom Prom	Ser Real Time	ial No.		s check	Internal standards Printer	
Eeprom		ial No. e clock	Temp	s check		
Eeprom Prom	Real Time Auto calib	ial No. e clock	Temp Motor Range	op:51	Printer	

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No:	68-012337	Ac	ccuracy check date/time: 06/13/2015	03:17
	bration check: bration check:		Interferent detected: Interferent detected:	
Concentration (%W/\	V) Solution lot#:	;	Soln temp C: Simulator Ser No.	
Invalid sample Number of people tes Operator	sted since last periodic dete	ermination of accu	uracy:	
		Diagnostics che	eck	
Ram Check sum 10489	Serial No. Real Time clock Auto calibration Download date/time: 0 d imported (Y/N) N	Temp Motor Range 06/15/2015 09:5 Source Ex		
Instrument serial No:	68-012337	Ac	ccuracy check date/time: 06/15/2015	09:55
	bration check: bration check: V) Solution lot#:		Interferent detected: Interferent detected: Soln temp C: Simulator Ser No.	
Number of people tes	sted since last periodic dete	ermination of accu	uracy:	
		Diagnostics che	eck	
Ram Pass , Check sum 23748	Serial No. Real Time clock Auto calibration Download date/time: 0 d imported (Y/N) N	Temp Fail Motor Range 06/15/2015 10:0 Source Di		

This is a true, complete, and accurate report of the periodic determination of accuracy and diagnostic data from the specified Intoxilyzer 5000EN instrument obtained from the CMI COBRA database.

5

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: 68-012337 Accuracy check date/time: 06/15/2015 09:57

First calibration check: Interferent detected: Second calibration check: Interferent detected:

Concentration (%W/V) Solution lot#: Soln temp C:

Simulator Ser No.

Number of people tested since last periodic determination of accuracy:

Operator

Diagnostics check

Internal standards Pass Serial No. Pass Temp Pass Eeprom Pass Prom Printer Pass **Pass** Real Time clock Pass Motor Pass Ram Pass Digital voltmeter 3231 Auto calibration Pass Range Pass Check sum 34097 Record ID 81J0 Download date/time: 06/15/2015 10:05

Record imported (Y/N) N Source Diagnostic

Instrument serial No: 68-012337 Accuracy check date/time: 06/15/2015 09:58

First calibration check:Not performed Interferent detected: No Second calibration check:Not performed Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C:

Simulator Ser No.

6

Number of people tested since last periodic determination of accuracy: 033

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 53325 Download date/time: 06/15/2015 10:05 Record ID 81G7

Record imported (Y/N) N Source Calibration

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: 68-012337 Accuracy check date/time: 06/15/2015 10:00

First calibration check: .112 Interferent detected: No Second calibration check: .116 Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C: 34.02

Simulator Ser No. DR3707

Number of people tested since last periodic determination of accuracy: 000

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 73130 Download date/time: 06/15/2015 10:05 Record ID 81H8

Record imported (Y/N) N Source Calibration

Instrument serial No: 68-012337 Accuracy check date/time: 06/18/2015 08:59

First calibration check: Interferent detected: Second calibration check: Interferent detected:

Concentration (%W/V) Solution lot#: Soln temp C:

Simulator Ser No.

7

Number of people tested since last periodic determination of accuracy:

Operator

Diagnostics check

Serial No. Pass Internal standards Pass Eeprom Temp Pass Pass Prom **Pass** Real Time clock Pass Motor Pass Printer Pass Ram **Pass** Range Pass Digital voltmeter 3443 Auto calibration Pass Check sum 51485 Record ID 81T3 Download date/time: 06/25/2015 09:42

Record imported (Y/N) N Source Diagnostic

This is a true, complete, and accurate report of the periodic determination of accuracy and diagnostic data from the specified Intoxilyzer 5000EN instrument obtained from the CMI COBRA database.

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: 68-012337 Accuracy check date/time: 06/18/2015 09:04

First calibration check: Not performed Interferent detected: No Second calibration check: Not performed Interferent detected: No

Concentration (%W/V) Solution lot#: Soln temp C:

Simulator Ser No.

Number of people tested since last periodic determination of accuracy: 004

Operator

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 39021 Download date/time: 06/25/2015 09:42 Record ID 81N0

Record imported (Y/N) N Source Calibration

Instrument serial No: 68-012337 Accuracy check date/time: 06/18/2015 09:08

First calibration check: 109 Interferent detected: No Second calibration check: 113 Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C: 34.03

Simulator Ser No. DR3707

8

Number of people tested since last periodic determination of accuracy: 004

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 95210 Download date/time: 06/25/2015 09:42 Record ID 8105

Record imported (Y/N) N Source Calibration

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: 68-012337 Accuracy check date/time: 06/18/2015 09:10

First calibration check: .115 Interferent detected: No Second calibration check: .116 Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C: 34.02

Simulator Ser No. DR3707

LAST CAL CHECK - COLD SIM HOSE

Number of people tested since last periodic determination of accuracy: 000

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 24422 Download date/time: 06/25/2015 09:42 Record ID 81P3

Record imported (Y/N) N Source Calibration

Instrument serial No: 68-012337 Accuracy check date/time: 06/20/2015 05:21

First calibration check: Interferent detected: Second calibration check: Interferent detected:

Concentration (%W/V) Solution lot#: Soln temp C:

Simulator Ser No.

9

Invalid sample

Number of people tested since last periodic determination of accuracy:

Operator

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 33584 Download date/time: 06/25/2015 09:42 Record ID 81V5

Record imported (Y/N) N Source Exception

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

IIIIOXII	yzer soootin Periodic	Determination of A	ccuracy / Maintenance Reco	Jiu
Instrument serial N	No: 68-012337	Accu	racy check date/time: 06/25/2015	5 09:19
First c	alibration check:	In	terferent detected:	
Second of	alibration check:	Ir	terferent detected:	
Concentration (%\	N/V) Solution		Soln temp C: mulator Ser No.	
Number of people Operator	tested since last periodic	determination of accura	cy:	
		Diagnostics chec	k	
Eeprom Pass Prom Pass Ram Pass	Serial No. Pass Real Time clock Pass Auto calibration Pass	Temp Pass Motor Pass Range Pass	Internal standards Pass Printer Pass Digital voltmeter 3427	
Check sum 39778	201111000000000000000000000000000000000			
Rec	cord imported (Y/N) N	Source Diag	gnostic	
Instrument serial N	No: 68-012337	Accu	racy check date/time: 06/25/2015	5 09:21
First o	alibration check: Not perfo	ormed In	terferent detected: No	
Second of	calibration check:Not perfo	ormed In	terferent detected: No	
Concentration (%\	N/V) Solution		Soln temp C: mulator Ser No.	
Number of people Operator	tested since last periodic	determination of accura	cy: 016	
		Diagnostics chec	<u></u> k	
Eeprom Prom Ram	Serial No. Real Time clock Auto calibration	Temp Motor Range	Internal standards Printer Digital voltmeter	
Check sum 55857		-	Record ID 81Q8	

This is a true, complete, and accurate report of the periodic determination of accuracy and diagnostic data from the specified Intoxilyzer 5000EN instrument obtained from the CMI COBRA database.

10

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: 68-012337 Accuracy check date/time: 06/25/2015 09:26

First calibration check: 109 Interferent detected: No Second calibration check: 113 Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C: 34.03

Simulator Ser No. DR3707

Number of people tested since last periodic determination of accuracy: 016

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 36224 Download date/time: 06/25/2015 09:42 Record ID 81R5

Record imported (Y/N) N Source Calibration

Instrument serial No: 68-012337 Accuracy check date/time: 06/25/2015 09:29

First calibration check: .114 Interferent detected: No Second calibration check: .115 Interferent detected: No

Concentration (%W/V) .120 Solution lot#: IRS121014 Soln temp C: 34.05

Simulator Ser No. DR3707

11

LAST CAL CHECK - COLD SIM HOSE

Number of people tested since last periodic determination of accuracy: 000

Operator G GOSSAGE 7885

Diagnostics check

EepromSerial No.TempInternal standardsPromReal Time clockMotorPrinterRamAuto calibrationRangeDigital voltmeter

Check sum 52082 Download date/time: 06/25/2015 09:42 Record ID 81S0

Record imported (Y/N) N Source Calibration

This is a true, complete, and accurate report of the periodic determination of accuracy and diagnostic data from the specified Intoxilyzer 5000EN instrument obtained from the CMI COBRA database.

1400 Canal Ave Long Beach, Ca 90813

PH: (562) 570-7680 Fax: (562) 570-7704

Intoxilyzer 5000EN Periodic Determination of Accuracy / Maintenance Record

Instrument serial No: Accuracy check date/time: 06/30/2015 68-012337 18:26 First calibration check: Interferent detected: Second calibration check: Interferent detected: Concentration (%W/V) Solution lot#: Soln temp C: Simulator Ser No. Number of people tested since last periodic determination of accuracy: Operator Diagnostics check Serial No. Pass Internal standards Pass Eeprom Pass Temp Pass Prom Printer Pass Pass Real Time clock Pass Motor Pass Ram Pass Auto calibration Pass Range Pass Digital voltmeter 3332

Check sum 90950 Download date/time: 07/02/2015 09:02

5 09:02 Record ID 82H9

Record imported (Y/N) N Source Diagnostic